OKATEST

ZE/OA48C

ZE/OA96C

Test for detection of Okadaic Acid-toxins group

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info@zeulab.com • www.zeulab.com

SCOPE

This protocol specifies a method for the quantitative determination of Okadaic Acid (OA) and other carboxylic toxins of the OA group including DTX1, DTX2 and DTX3 by a colorimetric phosphatase inhibition assay. This method is applicable to shellfish species such as mussels, clams, oysters and scallops.

PRINCIPLE

Test based on the phosphatase activity inhibition by OA-toxins group, responsible for diarrhetic shellfish poisoning (DSP).

Phosphatase enzyme PP2A is able to hydrolyse a specific substrate, yielding a product that can be detected at 405 nm. Samples containing toxins from the okadaic acid group will inhibit the enzyme activity proportionally to the amount of toxin contained in the sample. The concentration of toxin in the sample can be calculated using a standard curve.

KIT CONTENTS

	48 Tests Kit	96 Tests Kit
Microtiter plate strips (8 wells per strip)	6	12
Vials of Phosphatase (Phosphatase)	2	4
Set of Okadaic Acid Standards (Okadaic acid 0.5, 0.8, 1.2, 1.8 and 2.8 nM)	1	1
Chromogenic Substrate (Chromogenic Substrate)	1	1
Phosphatase Dilution Buffer (Phosphatase Dilution Buffer)	1	1
Stock Buffer Solution (Stock Buffer Solution)	1	1
Stop Solution (Stop Solution)	1	1
Adhesive film	1	2
Kit instructions	1	1

ADDITIONAL MATERIAL AND REAGENTS NEEDED

- Micropipettes
- Blender (Ultraturax) or mortar and pestle
- Heater at $30^{\circ}\text{C} \pm 2^{\circ}\text{C}$ (i.e. FX Incubator, Ref ZE/FX, from ZEU-INMUNOTEC)
- Microplate reader (wavelength at 405 nm)
- Water bath for $76 \pm 2^{\circ}\text{C}$
- Methanol (analytical grade)
- NaOH 2.5 N made by titration, (NaOH of analytical grade)
- HCl 2.5 N made by titration, (HCl of analytical grade)
- Deionised water (grade 2, ISO3696)
- Graded 50 mL centrifuge tubes with screw caps
- Tube shaker

SOLUTIONS

- 1.- **Okadaic Acid Standards:** to make sure these solutions are homogeneous, it is very important to mix well using a vortex, before applying to the plate.
- 2.- **Chromogenic Substrate solution:** The solution contains stabilization resin. Make sure this resin is not added to the microwells. To assure that, it is recommended to transfer the volume needed into a transparent labware (i.e.: test tube or eppendorf) and take the solution from that container to add into the wells. *Note:* Do not use this solution if the absorbance of 90 μ L of this solution at 405 nm is over 0.6.
- 3.- **Phosphatase solution:** Add 2.0 mL of phosphatase dilution buffer (**Phosphatase Dilution Buffer**) to one of the phosphatase vials (**Phosphatase**) and dissolve by mixing gently for 1 hour $\pm$ 5 minutes at room temperature (22 ± 2 °C) to ensure that the enzyme is fully hydrated. **Do not use the tube shaker at any moment.** This solution must be stored under refrigeration if not in use immediately after preparation. Do not use the phosphatase solution for following days. Each enzyme vial contains enough volume for 24 wells. If more than one vial is used in the assay, dissolve each vial as described above, make a pool with the content of the vials and mix gently, by inversion, before use.
***Attention:** this reagent is blue and becomes brownish when dissolved. If brownish colour is noticed before hydration, discard this reagent as it could be damaged.
- 4.- **Buffer solution x1:** dilute the **Stock Buffer Solution** included in the kit by mixing 1 volume with 9 volumes of deionised water. Use buffer solution x1 only freshly made, and store under refrigeration if not in use immediately after preparation.
- 5.- **2.5 N NaOH:** weigh 100 g of NaOH and add 500 mL of water and dissolve. Transfer to a volumetric flask and add deionised water up to a final volume of 1000 mL.
- 6.- **2.5 N HCl:** add 205 mL of HCl (37 %) to 400 mL of deionised water already contained in a volumetric flask. Make the volume up to 1000 mL with deionised water.

SAMPLES EXTRACTION

The method described below includes a hydrolysis step to detect all toxins forms of okadaic acid (okadaic acid and dinophisistoxins).

- 1.- Clean the shell thoroughly using water
- 2.- Open the shellfish by cutting the adductor muscles.
- 3.- Wash inside the shell thoroughly to get rid of any dirt.
- 4.- Remove the tissue inside the shell by cutting all the muscles attached to the shell.
- 5.- Place the shellfish tissue in a filter paper for few minutes to remove water in excess.

It is recommended to use graded 50 mL centrifuge tubes with screw caps during the following steps of hydrolysis in order to prevent losses due to labware changes.

- 6.- Mash the shellfish tissue to obtain a representative sample and weigh 5 g. Add 25 mL of Methanol and homogenise the mixture for 2 minutes using a tube shaker.
- 7.- Centrifuge at 2000 g for 10 min at 4 °C. The supernatant (*methanolic extract*) is poured into a centrifuge tube.
- 8.- Take 640 μ L of *methanolic extract* and pour into another centrifuge tube.
- 9.- Add 100 μ L of 2.5 N NaOH.
- 10.- Seal and heat at 76 ± 2 °C for 40 minutes.

- 11.- Add 80 μ L of 2.5 N HCl (the sample does not need to be cooled down previously).
- 12.- Add up to 20 mL of **Buffer solution x1**.

TEST PROCEDURE

Warning:

The volume of some reagents used in this assay is small and special attention must be paid when added to the wells:


- Make sure the pipettes are calibrated before running the assay.
- Use pipettes according to the volumes to be dispensed. Use pipettes with a maximum pipette volume of 100 or 200 μ L.
- Be sure that the incubator's temperature is stabilized before use.

It is recommended to run samples and standards in duplicate.

- 1.- Add 50 μ L of samples or standards.
- 2.- Add 70 μ L of the Phosphatase Solution to each well. Mix well by gentle tapping on the side of the plate.
- 3.- Cover the plate with the adhesive film provided and incubate for 20 ± 0.5 minutes at 30 ± 2 °C.
- 4.- Remove the adhesive film and add 90 μ L of Chromogenic Substrate to each well. Mix well by gently tapping on the side of the plate.
- 5.- Cover the plate with the adhesive film and incubate 30 ± 0.5 minutes at 30 ± 2 °C.
- 6.- Remove the adhesive film and add 70 μ L of Stop Solution to each well.
- 7.- Read absorbance of samples and standards at 405 nm.

GRAPHIC REPRESENTATION AND CALCULATIONS OF RESULTS

- 1.- Obtain a standard curve by plotting the absorbance values in a linear *y axis* and the concentration of okadaic acid in a logarithmic *x axis* and use a logarithmic fitting as shown in the graphic next page. R^2 has to be greater than or equal to 0.96.

- 2.- The OA concentration contained in the sample (C_s) is calculated by interpolation into the calibration curve or using the following equation:

$$x = \text{EXP}(y - b/a)$$

Where x is the OA concentration in the sample (C_s) and y the absorbance of the sample.

Note: An Excel worksheet to calculate results is available upon request.

- 3.- Calculate the diarrhetic shellfish toxins concentration in tissue (C_t) as follows:

$$C_t (\mu\text{g/kg}) = \frac{(C_s (\text{nM}) \times \text{FD} \times \text{MW} (\text{g/mol}) \times V_e (\text{L}))}{M_t (\text{g})}$$

C_t: toxins concentration in tissue, expressed as equivalents of OA; **C_s:** toxins concentration in sample; **FD:** Methanolic extract dilution factor (i.e. 640 $\mu\text{L}/20 \text{ mL} \rightarrow \times 31.25$); **MW:** Okadaic acid molecular weight = 805; **V_e:** Methanolic extract volume (0.025L); **M_t:** Tissue weight (5g).

Example: for OA concentration of 1.5 nM: $1.5 \text{ nM} \times 31.25 \times 805 \text{ g/mol} \times 0.025 \text{ L} / 5 \text{ g} = 189 \text{ } \mu\text{g OA eq/kg}$.

For samples with OA concentration falling outside the working range ($< 0.5 \text{ nM}$ or $> 2.8 \text{ nM}$), results will be reported as $< 0.5 \text{ nM}$ (or $< 63 \text{ } \mu\text{g/Kg}$) or $> 2.8 \text{ nM}$ (or $> 352 \text{ } \mu\text{g/kg}$), respectively.

STABILITY AND STORAGE

The kit contents must be stored at $4 - 12 \text{ }^{\circ}\text{C}$ and protected from light. This kit has a shelf life of 4 - 6 months when stored under optimal conditions. See the expiry date on the kit package.

SAFETY

Safety clothing should be worn and skin contact with the reagents avoided. Do not ingest. A SAFETY DATA SHEET is available from your local distributor on request.

***Warning:** Okadaic Acid is toxic. Gloves, mask and other protective clothing must be worn when handling okadaic acid solutions.

REFERENCES

1. A Fluorescent Microplate Assay for Diarrheic Shellfish Toxins. Vieytes M. R. et al. Analytical Biochemistry 248, 258-264 (1997).
2. Protein phosphatase inhibition assay adapted for determination of total DSP in contaminated mussels. Mounfort D. O. et al. Toxicon 39, 383-390 (2001).
3. Inter-laboratory validation of the fluorescent protein phosphatase inhibition assay to determine diarrheic shellfish toxins: intercomparison with liquid chromatography and mouse bioassay. Gonzalez J. C., et al. Analytica Chimica Acta 466, 233-246 (2002).

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
- 1.- A Fluorescent Microplate Assay for Diarrheic Shellfish Toxins. Vieytes M. R. et al. Analytical Biochemistry 248, 258-264 (1997).
- 2.- Protein phosphatase inhibition assay adapted for determination of total DSP in contaminated mussels. Mounfort D. O. et al. Toxicon 39, 383-390 (2001).
- 3.- Inter-laboratory validation of the fluorescent protein phosphatase inhibition assay to determine diarrheic shellfish toxins: intercomparison with liquid chromatography and mouse bioassay. González J. C., et al. Analytica Chimica Acta 466, 233-246 (2002).

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